

Experiment 5

Identification of a Metal Carbonate

Pre-Lab Assignment

Before coming to lab:



- Read the lab thoroughly.
- Answer the pre-lab questions that appear at the end of this lab exercise.

Purpose

The chemical formula of an unknown metal carbonate or bicarbonate will be identified from a given list by its molar mass. The molar mass will be determined via stoichiometry after adding acid to produce and lose carbon dioxide gas.

Background

Carbonate (CO_3^{2-}) and bicarbonate (HCO_3^-) are polyatomic ions that pair with metals to form ionic compounds. When a metal carbonate or bicarbonate is reacted with hydrochloric acid, a water-soluble ionic compound, water, and carbon dioxide form. This is the "acid test" often used by geologists to identify carbonates in rocks and follows the balanced reaction in Eqn. 1 for a metal carbonate and Eqn. 2 for a metal bicarbonate. In these equations, M = Na or K.



While carbon dioxide is fairly soluble in water under high pressures, as evidenced by canned soft drinks, it escapes at normal atmospheric pressures. This explains why when left open, containers of soda go flat rather quickly. Carbon dioxide is slightly more dense than atmospheric air so will sink.

Stoichiometry is the study of amounts. In chemistry, stoichiometry allows us to calculate the amount of a particular chemical required for a reaction and how it relates to the other reactants and products. For this experiment, the metal carbonate sample will be fully reacted with hydrochloric acid to produce the products seen in Eqn. 1 and 2. By the Law of Conservation of Mass, the mass of the reactants before reaction should equal the mass of products after. However, since $\text{CO}_2(\text{g})$ is a gas, the container will be left open and allowed to escape so that the amount of $\text{CO}_2(\text{g})$ produced can be measured by subtraction. Using the stoichiometric relationships in Eqn. 1 and 2 can calculate the amount of moles of the original metal carbonate that reacted.

The molar mass of a substance is usually determined from the atomic masses on the periodic table. However, since the exact formula of the metal carbonate is unknown, the molar mass will be experimentally determined as the grams of metal carbonate reacted/moles of metal carbonate reacted and then compared to the options listed on Table 1 to identify its formula.

Table 1: Possible Metal Carbonate Unknowns

Chemical Name	Chemical Formula	Molar Mass
sodium carbonate	Na ₂ CO ₃	105.988 g/mol
sodium carbonate · ½ hydrate	Na ₂ CO ₃ · ½ H ₂ O	114.996 g/mol
sodium bicarbonate	NaHCO ₃	84.066 g/mol
potassium carbonate	K ₂ CO ₃	138.206 g/mol
potassium bicarbonate	KHCO ₃	100.115 g/mol
potassium carbonate · 1 ½ hydrate	K ₂ CO ₃ · 1 ½ H ₂ O	165.229 g/mol

Example Problem: Determining the Molar Mass of an Unknown

An empty 125 mL Erlenmeyer flask weighs 85.0155 g. An unknown metal carbonate is added to the flask and weighed again at 86.0255 g. It and 15 mL of 6.0 M HCl in a styrofoam weigh 101.1115 g. The HCl solution is added to the flask, allowed to react, and the CO₂(g) escaped. The resulting mixture after the reaction weighed 100.7250 g. Find the molar mass of the unknown and identify it on Table 1.

Step 1: Find the mass of metal carbonate

$$86.0255 \text{ g} - 85.0155 \text{ g} = 1.0100 \text{ g}$$

Step 2: Find the mass of CO₂(g) that escaped

$$101.1115 \text{ g} - 100.7250 \text{ g} = 0.3865 \text{ g}$$

Step 3: Find the moles of CO₂(g)

$$0.3865 \text{ g CO}_2 \times \frac{1 \text{ mol CO}_2}{44.0095 \text{ g CO}_2} = 0.008782 \text{ moles CO}_2$$

Step 4: Find the moles of metal carbonate

$$0.008782 \text{ moles CO}_2 \times \frac{1 \text{ mol metal carbonate}}{1 \text{ mol CO}_2} = 0.008782 \text{ moles metal carbonate}$$

Step 5: Find the molar mass of the metal carbonate

$$\frac{1.0100 \text{ g}}{0.008782 \text{ moles}} = 115.1 \frac{\text{g}}{\text{mol}} \text{ sodium carbonate} \cdot \frac{1}{2} \text{ hydrate}$$

Procedure

1. Weigh a clean, dry 125 mL Erlenmeyer flask. Record this mass in your data sheet.
2. Add 1.0-2.0 g of your metal carbonate to the Erlenmeyer flask and reweigh. Record this mass in your data sheet.
3. Collect a clean, dry Styrofoam cup. Measure approximately 15 mL of 6.0 M HCl(aq) and pour it carefully into the cup, being sure to not spill.
4. Place the Erlenmeyer flask from Step 2 on the balance. Carefully set the cup from Step 3 on the top. You may need to open the top door on the balance. Do not let the cup or flask rest on any of the balance's doors. Record the mass in your data sheet for as many digits remain stable (the mass will continue to decrease due to the solution's evaporation).
5. At your bench, slowly pour the HCl into the Erlenmeyer flask. Do not spill.
6. After the bubbling subsides, swirl the flask and gently blow into the flask to push the CO₂(g) out the top. Do not inhale. Continue to do so until all the bubbling has stopped.
7. Place the Erlenmeyer flask and contents on the balance. Carefully set the empty cup on top and record the weight in your data sheet.
8. Dispose of the neutralized solution in the sink with running water. Wash and dry all equipment used.
9. Repeat Steps 1-8 for a second trial. Your two trials should have an average molar mass that matches with a value shown in Table 1.

Do not allow your flask's contents to foam over!

HCl releases corrosive fumes that should not be inhaled.

Experiment 5—Data Sheet

Name: _____

Unknown #: _____

	Trial One	Trial Two
1. Mass Erlenmeyer flask (g)	_____	_____
2. Mass Erlenmeyer flask and unknown (g)	_____	_____
3. Mass flask, unknown, cup, and HCl before reaction (g)	_____	_____
4. Mass flask, unknown, cup, and HCl after reaction (g)	_____	_____
5. Mass unknown (g) <i>show calculation:</i>	_____	_____
6. Mass CO ₂ (g) lost (g) <i>show calculation:</i>	_____	_____
7. Moles CO ₂ (g) (mols) <i>show calculation:</i>	_____	_____
8. Moles unknown (mols) <i>show calculation:</i>	_____	_____
9. Molar Mass unknown (g/mol) <i>show calculation:</i>	_____	_____
10. Average Molar Mass of unknown (g/mol)	_____	_____
11. Identity of unknown (from Table 1):	_____	_____

Experiment 5—Post-Lab Assignment

1. Calculate the percent error between your average molar mass and the tabulated value from Table 1.

2. Give three experimental reasons why your average molar mass may differ from the tabulated value.

3. A student performed the identical experiment as described to determine the identity of an unknown metal carbonate. They completed every step perfectly *except* they did not blow out the dense $\text{CO}_2(\text{g})$ after the reaction. Would this change the recorded mass of the flask, sample, cup, and HCl after reaction? If so, would the new value be higher or lower? Explain.

4. How would the error described in Question 3 change (too high, too low, no effect) the calculated molar mass of the unknown? Explain.

5. A student performed the identical experiment as described to determine the identity of an unknown metal carbonate. They completed every step perfectly *except* they had a small piece of glass in the Erlenmeyer flask that remained from the beginning and throughout the entire experiment. Would this error make the calculated molar mass of the unknown inaccurate? If so, would the incorrect value be higher or lower than the correct one? Explain.

Experiment 5—Pre-Lab Assignment

Name: _____

For all calculations, show all work and draw a box around the final answers.

1. Why are you instructed to use no less than 1.0 g and no more than 2.0 g of sample?

2. A student obtained the following data from this experiment:

Mass flask (g)	75.0235 g
Mass flask + unknown (g)	76.8912 g
Mass flask, unknown, cup, and HCl before reaction	101.5894 g
Mass flask, unknown, cup, and HCl after reaction	102.3649 g

a. Calculate the molar mass in g/mol of the student's unknown.

b. According to Table 1, which metal carbonate is the unknown?

c. Write the balanced chemical equation for the student's unknown and the reaction performed in this experiment.